

QwikCheck™ QC EQUINE

Sperm Quality Analyzer

QC an EQUINE AI sample
prior to insemination in two minutes!

Fast, simple and accurate!



- Target market: Horse owners/VETs who inseminate mares and want to QC semen samples
- Tests EXTENDED, COOLED and FROZEN equine semen
- Print results to a label maker
- Re-usable, washable testing capillary

SEMEN parameters Reported:

- Sperm Concentration (Not reported for Frozen)
- Motility
- Progressive Motility
- Motile Sperm Concentration
- Prog. Motile Sperm Concentration
- Velocity



"Remember, it all started with a Sperm!"